



S I N G U L A R  
G E N O M I C S



# G4X VIEWER



# G4X FILE OUTPUTS

# Understanding G4X File Outputs Names

## Possible Tissue Locations

4-lane flow cell

A1	B1	C1	D1	E1	F1	G1	H1
A2	B2	C2	D2	E2	F2	G2	H2
A3	B3	C3	D3	E3	F3	G3	H3
A4	B4	C4	D4	E4	F4	G4	H4

2-lane flow cell

A1	B1	C1	D1	E1
A2	B2	C2	D2	E2

Name	Status	Date	Type	Size
g4-012-025-FC2-L001_C01.bin	✓	1/24/2025 11:34 AM	BIN File	944,979 KB
g4-012-025-FC2-L001_C01.ome	✓	1/24/2025 12:43 PM	TIFF File	33,591,097 KB
g4-012-025-FC2-L001_C01.tar	✓	1/24/2025 11:34 AM	TAR File	3,147,850 KB
g4-012-025-FC2-L001_C01_HE.ome	✓	1/24/2025 11:34 AM	TIFF File	4,812,201 KB
g4-012-025-FC2-L001_C01_nuclear.ome	✓	1/24/2025 11:35 AM	TIFF File	4,540,779 KB

Instrument name/  
number


Run number

Flow cell number  
\* The nest can take  
FC1, FC2, FC3, FC4  
depending on where  
the sample is placed

Lane number

Tissue location on the  
lane

# G4X File outputs




Name	Type
C03.ome	TIFF File
C03	Compressed Archive ...
C03_HE.ome	TIFF File
C03_nuclear.ome	TIFF File
C03_run_metadata	JSON Source File
colon_cancer_rep1.bin	BIN File

File output extension	Data type
<code><sample_1>.bin</code>	Cell segmentation mask
<code><sample_1>.ome</code>	Proteomics data
<code><sample_1>.tar</code> folder	Transcriptomics data
<code><sample_1>_HE.ome</code>	fH&E images
<code><sample_1>_nuclear.ome</code>	Fluorescent version of the fH&E
<code><sample_1>_metadata.json</code>	General Details File

# G4X VIEWER

# G4X HOMEPAGE



G4X Viewer Version

Multifile Upload  
Upload each individual  
file

Single file Upload  
requires a compressed .tar  
folder

html-based G4X viewer : <https://g4x-viewer.singulargenomics.com/>

# UPLOADING FILES

The screenshot shows the G4X Viewer web application interface. The main content area displays the text "Please upload an image file to view." On the right side, there is a sidebar with the G4X Viewer logo and version number "2.6.1". The sidebar contains several upload sections:

- Image File Name:** A text input field with an upload icon (cloud with arrow) to its right.
- UPLOAD IMAGE FILE:** A button with a dashed border.
- General Details File:** A text input field.
- UPLOAD GENERAL DETAILS:** A button.
- Transcript File Name:** A text input field.
- UPLOAD POINTS FILE:** A button.
- Cell Mask Filename:** A text input field.
- UPLOAD CELL MASK:** A button.

At the bottom of the sidebar, there are social media icons for a globe, LinkedIn, X, GitHub, and Email.

Proteins images file\*  
Ome.tiff

Metadata JSON file  
.json

Transcripts Images Folder  
.tar



Cell Segmentation Mask  
.bin



# PROTEOMICS

# PROTEIN CHANNEL SETTINGS

 **G4X Viewer**



Mouse Pos: [1716, 7968] Layer: 5/6 Share

 **G4X Viewer**  
2.6.1

**SOURCE FILES**

MULTI-FILE  SINGLE-FILE

Image File Name  
C03.ome.tiff

**UPLOAD IMAGE FILE**

General Details File

**UPLOAD GENERAL DETAILS**

Transcript File Name

**UPLOAD POINTS FILE**

Cell Mask Filename

**UPLOAD CELL MASK**

**VIEW SETTINGS**

**PROTEIN CHANNEL SETTINGS**


**BRIGHTFIELD IMAGES SETTINGS**

# PROTEIN CHANNEL SETTINGS

Choose a Proteomics Color Map to visualize the proteins selected in a heat map to show protein density


Overview window at the bottom left – especially useful when visualizing a zoomed in anatomical region on the tissue

You can export and import your protein color and contrast settings to reuse or share for consistency across datasets or publications.



Overview window

# LENS FUNCTION



**G4X Viewer**  
2.6.2

Overview  
 Lens

Highlighted Channel:  
PanCK

Channel Controls

EXPORT IMPORT


nuclear  
274.0  
35 4244

PanCK  
87.0  
1 4578


CD3  
145.0  
1 2010

Can be turned on for a specific marker of interest to visualize individually on the tissue

Greys out all proteins except for the one highlighted



# PROTEOMICS CHANNEL SETTINGS



**G4X Viewer** 2.6.2

EXPORT IMPORT

nuclear 76.0 35 4244

PanCK 12.0 1 4578

αSMA 4.0 1 6531

CD31 3.0 1 7091

+ ADD CHANNEL

Change Color of individual markers

Pixel intensity of individual channel


Adjust contrast

Turn on and off individual channels

Visualize up to 6 Channels

# GENERAL DETAILS FILE

# GENERAL DETAILS FILE




**Run Metadata**

Machine:	g4-007
Run Id:	127
Platform:	g4x
Fc:	FC3
Lane:	L003
Time Of Creation:	2025-03-18 23:51:10
Transcript Panel:	["iss-t-0242", "iss-t-0243", "iss-t-0244", "iss-t-0245", "iss-t-0246", "iss-t-0247"]
Protein Panel:	["iss-p-0024", "iss-p-0025", "iss-p-0026", "iss-p-0027", "iss-p-0028"]
Software:	caretta
Software Version:	2.11.1

**Core Metrics**

Sample Id:	C03
Software:	caretta
Software Version:	2.11.1
Total Tissue Area Mm <sup>2</sup> :	16.958955078125
Number Cells:	89820
Median Cell Area Um <sup>2</sup> :	107.51953125
Median Transcripts Per Cell:	72
Median Unique Genes Per Cell:	25
Pct Empty Cells:	0.3017145401914941
Ncp Fdr Pct:	0.0465591752731603
Ncs Fdr Pct:	0.010399999687453
Pct Transcripts In Cells:	81.12692293012093
Total Transcripts:	11091445
Total Ncp:	87
Total Ncs:	44
Median Transcripts Per 100um <sup>2</sup> :	50
Pct Q20:	91.53835036699172



**G4X Viewer** 2.6.1

**SOURCE FILES**

MULTI-FILE SINGLE-FILE

Image File Name  
C03.ome.tiff

UPLOAD IMAGE FILE

General Details File  
C03\_run\_metadata.json

UPLOAD GENERAL DETAILS

Transcript File Name

UPLOAD POINTS FILE

Cell Mask Filename

UPLOAD CELL MASK

**VIEW SETTINGS**

**PROTEIN CHANNEL SETTINGS**

**BRIGHTFIELD IMAGES SETTINGS**

**TRANSCRIPT LAYER SETTINGS**

**SEGMENTATION LAYER SETTINGS**

Once the metadata.json file is uploaded, the info icon becomes active, showing assay and sample details

Upload the metadata.json file

# TRANSCRIPTS


# TRANSCRIPTS: UPLOAD .TAR FOLDER

The screenshot shows the G4X Viewer interface. On the left, a large visualization displays a spatial transcriptomics dataset with colored dots representing individual cells. A smaller inset shows a zoomed-in view of a specific region. The right sidebar contains the 'G4X Viewer 2.6.1' control panel. It features a 'SOURCE FILES' section with two tabs: 'MULTI-FILE' (selected) and 'SINGLE-FILE'. Below the tabs are three upload sections: 'Image File Name' (C03.ome.tiff) with an 'UPLOAD IMAGE FILE' button; 'General Details File' (C03\_run\_metadata.json) with an 'UPLOAD GENERAL DETAILS' button; and 'Transcript File Name' (C03.tar) with an 'UPLOAD POINTS FILE' button. A pink dashed box highlights the 'Transcript File Name' section, and a pink arrow points from the text 'Upload .tar compressed folder' to the 'UPLOAD POINTS FILE' button.

Name	Type
C03_nuclear.ome	TIFF File
C03_HE.ome	TIFF File
C03.ome	TIFF File
C03_run_metadata	JSON Source File
C03	Compressed Archive Folder
colon_cancer_rep1.bin	BIN File

Upload .tar compressed folder

# VIEW SETTINGS



**G4X Viewer**  
2.6.2

**VIEW SETTINGS**

**Global Selection**

T

Z

**Zoom Control**

Scale

**Layers Toggles**

Transcript Layer

Show Tile Boundaries

Cell Masks Layer

**PROTEIN CHANNEL SETTINGS**

In version 2.6.2 T and Z are not functional but will be updated in later versions

Can set the scale manually

Toggle the transcript layer on and off

Toggle the cell mask layer on and off

# TRANSCRIPT LAYER SETTINGS

The screenshot shows the G4X Viewer interface. On the left is a large visualization of a cell with numerous blue and orange points representing transcripts. On the right is a settings panel for the 'TRANSCRIPT LAYER SETTINGS'. The settings include:

- BRIGHTFIELD IMAGES SETTINGS** (collapsed)
- TRANSCRIPT LAYER SETTINGS** (expanded)
- Advanced Options**
  - Enable layers controls
  - Percent of transcripts shown: 100%, 20%, 4%, 0.8%, 0.16% (slider)
- Point Size**
  - 1.5 (slider)
- Point Filters**
  - Enable Filter
  - Show Discarded
- Gene List**

GENE NAME	COLOR
ABCG2	Blue
ACE2	Green
ACKR1	Red
ACTA2	Brown
ACTG2	Olive
ADIPOQ	Yellow
AGR3	Teal
- SEGMENTATION LAYER SETTINGS** (collapsed)

At the bottom of the visualization, there is a mouse position indicator: Mouse Pos: [9226, 11402] Layer: 5/6 Shape: 1, 1, 16, 1200, 952. A scale bar indicates 348.1 μm.

Enable Layers allows you to choose the percent of transcripts you want to view


Change the point size visible genes

By default, you are viewing all genes. By enabling filter, you can select only the genes you want to view

List of genes and the color they are being displayed


\*Note the default is 0.16% of all transcripts. When you zoom in the percent of transcripts visible increases. The max zoom shows 100% of transcripts

# TRANSCRIPTS: ENABLE LAYERS CONTROL




When Enabling layers controls a warning message will pop up to let you know that viewing 100% of transcripts can crash the application

You can confirm to continue



Now you can choose the percent of genes you want to view. Note that it does take time to render

# TRANSCRIPTS: ENABLE FILTER



Enable filter shows only selected genes

Search for specific genes to visualize

Select to view only genes you have selected





fH&E

# BRIGHTFIELD IMAGES SETTINGS

The screenshot displays the G4X Viewer interface. On the left, a large brightfield image of a tissue section is shown with overlaid green and blue channels. A smaller thumbnail of the same image is visible in the bottom-left corner. At the bottom of the image, a scale bar indicates 355.7 μm, and the mouse position is [14294, 2950]. The layer information is 1/6, and the shape is 1, 1, 16, 19200, 15232. On the right, the settings sidebar is open, showing the G4X Viewer logo and version 2.6.2. The sidebar is divided into several sections: SOURCE FILES, which includes upload buttons for MULTI-FILE and SINGLE-FILE, and input fields for Image File Name (C03.ome.tiff), General Details File, Transcript File Name, and Cell Mask Filename, each with an associated UPLOAD button; VIEW SETTINGS; PROTEIN CHANNEL SETTINGS; BRIGHTFIELD IMAGES SETTINGS (highlighted with a pink oval); TRANSCRIPT LAYER SETTINGS; and SEGMENTATION LAYER SETTINGS. At the bottom of the sidebar, there are social media icons for globe, LinkedIn, X, GitHub, and email.

# BRIGHTFIELD IMAGES SETTINGS: UPLOAD fH&E



The opacity of the fH&E image can be adjusted from here to make it easier to switch back and forth between proteins and fH&E

Upload fH&E

\*fH&E image can be visualized on top of the proteins and the transcripts data

# CELL SEGMENTATION


# CELL LAYER SETTINGS: UPLOAD.BIN FILE

The screenshot displays the G4X Viewer interface. On the left, a large image shows a cell segmentation mask with various colored regions (blue, green, orange, purple) on a black background. The right sidebar contains the settings panel, which is divided into several sections: SOURCE FILES, VIEW SETTINGS, PROTEIN CHANNEL SETTINGS, BRIGHTFIELD IMAGES SETTINGS, TRANSCRIPT LAYER SETTINGS, and SEGMENTATION LAYER SETTINGS. The SEGMENTATION LAYER SETTINGS section is highlighted with a pink dashed box. Within this section, the 'Cell Mask Filename' field is set to 'colon\_cancer\_rep1.bin', and the 'UPLOAD CELL MASK' button is also highlighted with a pink dashed box. A pink arrow points from the text 'Upload cell segmentation mask - .bin file' to the 'UPLOAD CELL MASK' button. Another pink arrow points from the text 'Segmentation Layer Settings' to the 'SEGMENTATION LAYER SETTINGS' section header.

Upload cell segmentation mask - .bin file

Segmentation Layer Settings

# CELL LAYER SETTINGS



Toggle on and off cell segmentation

Can adjust the opacity of the cell mask




Image zoomed in to see cell mask

# CELL LAYER SETTINGS: ENABLE FILTER

G4X Viewer 2.6.1

VIEW SETTINGS

PROTEIN CHANNEL SETTINGS

BRIGHTFIELD IMAGES SETTINGS

TRANSCRIPT LAYER SETTINGS

SEGMENTATION LAYER SETTINGS

Cell Fill

Show cell fill

87

Cell Filters

Enable Filter

Show Discarded

Search...

<input checked="" type="checkbox"/>	CLUSTER ID	COLOR
<input checked="" type="checkbox"/>	3.0	Orange
<input checked="" type="checkbox"/>	6.0	Brown

1-2 of 2

Show active filters only

270.5 μm

Mouse Pos: [13106, 10352] Layer: 5/6 Shape: 1, 1, 16, 1200, 952

Enable filters to only view selected clusters

Search cluster IDs

Only view selected clusters

# CELL LAYER SETTINGS: UMAP

The image shows a software interface for G4X Viewer (version 2.6.2). The main view is a cell layer visualization with a dark background and colorful cell clusters. An interactive UMAP plot is overlaid on the left, showing a scatter plot of cell clusters in green. The UMAP plot has a grid and axes ranging from -15 to 20 on the x-axis and -20 to 15 on the y-axis. The UMAP plot window includes controls for Point Size (set to 1) and Subsampling Step (set to 1). Below the plot are input fields for X Start, X End, Y Start, and Y End, along with CLEAR and CONFIRM buttons. A scale bar at the bottom right indicates 454.1 μm. The mouse position is shown as [2065, 579] Layer: 5/6 Shape: 1, 1, 16, 1200, 952. On the right, a settings panel for G4X Viewer is visible. It includes a Cell Fill section with a 'Show cell fill' toggle and a slider set to 56.9. The Cell Filters section has 'Enable Filter' and 'Show Discarded' toggles. Below is a search bar and a table of clusters. The table has columns for CLUSTER ID and COLOR. The clusters listed are 3.0 (orange), 4.0 (green), 5.0 (purple), 6.0 (brown), 7.0 (grey), nan (yellow), and -1 (grey). There are 1-8 of 8 clusters shown. Below the table is a 'Show active filters only' checkbox. The Graph Filters section has a pink box around the UMAP filter, with FLOW CYTOMETRY below it. At the bottom of the settings panel are social media icons for globe, LinkedIn, X, GitHub, and email.

Umap

Point Size: 1 (Min: 1 | Max: 10) Subsampling Step: 1 (Min: 1 | Max: 20)

X Start: X End: Y Start: Y End: CLEAR CONFIRM

454.1 μm

Mouse Pos: [2065, 579] Layer: 5/6 Shape: 1, 1, 16, 1200, 952

G4X Viewer 2.6.2

Cell Fill: Show cell fill (56.9)

Cell Filters: Enable Filter Show Discarded

Search...

CLUSTER ID	COLOR
3.0	Orange
4.0	Green
5.0	Purple
6.0	Brown
7.0	Grey
nan	Yellow
-1	Grey

1-8 of 8



Show active filters only

Graph Filters: UMAP FLOW CYTOMETRY

In addition to visualizing individual cell clusters, you can utilize the UMAP selector tool to select a specific ROI on the interactive UMAP embedding and plot the spatial localization of selected cell.

This can be useful for exploring cluster- and sub-cluster level spatial patterning of cells in your tissue.


# CELL LAYER SETTINGS: UMAP FILTER



Spatial plotting of cells from selected ROI

- To use the selector tool, use the cursor to select a specific group of cells within the UMAP you want to spatially plot
- Once you have selected an ROI, hit confirm and close out of the window to update the viewer

# CELL LAYER SETTINGS: FLOW CYTOMETRY



The G4X Viewer 2.6.2 settings panel includes the following sections:

- Cell Fill:** A toggle for 'Show cell fill' is turned on, with a slider set to 56.9.
- Cell Filters:** A toggle for 'Enable Filter' is turned on, and 'Show Discarded' is turned off. A search bar is present above a table of filters.



<input type="checkbox"/>	CLUSTER ID	COLOR
<input type="checkbox"/>	3.0	Orange
<input type="checkbox"/>	4.0	Green
<input type="checkbox"/>	5.0	Purple
<input type="checkbox"/>	6.0	Brown
<input type="checkbox"/>	7.0	Grey
<input type="checkbox"/>	nan	Yellow
<input type="checkbox"/>	-1	Light Grey

Below the table is a 'Show active filters only' checkbox and a 'Graph Filters' section with two buttons: 'UMAP' and 'FLOW CYTOMETRY'. The 'FLOW CYTOMETRY' button is highlighted with a pink box. At the bottom are icons for globe, LinkedIn, X, refresh, and mail.

A Flow Cytometry ROI filter can also be applied to the data, enabling you to select cells based on protein expression levels. This filter can be applied independently or in conjunction with others like the UMAP ROI filter.

\*When you first open the Flow Cytometry window, it is in the linear scale

# CELL LAYER SETTINGS: FLOW CYTOMETRY FILTER




Confirm and close after making your selection to visualize cells

\*If you did not clear the UMAP this will only show cells selected within the UMAP filter

You can also select an ROI with the flow cytometry selection tool. Here, we are plotting all cells which are CD3 and CD8 positive.

# CELL LAYER SETTINGS: FLOW CYTOMETRY FILTER SETTINGS





Change graph settings

Change Color

Can change Axis Type

# CELL LAYER SETTINGS: FLOW CYTOMETRY FILTER LOG SCALE




Close and confirm to visualize cells selected  
\*If you did not clear the UMAP this will only show cells selected within the UMAP filter


Select cells of interest with selection tool.  
Here we chose CD3 and CD8 positive cells

# OTHER FEATURES


# OTHER FEATURES



When cell mask is enabled  
Hover over a cell to reveal detailed information



When the transcript layer is enabled  
Hover over a single transcript to reveal detailed information



The scale bar and screenshot matched current resolution of view

For more information find us online with the links in the bottom right corner

# G4X VIEWER



QUESTIONS?



YOU CAN CONTACT THE CARE TEAM AT  
[CARE@SINGULARGENOMICS.COM](mailto:CARE@SINGULARGENOMICS.COM)